Agitatori Magnetici


**MAGNETIC
MIXER**

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Principle of Operation

The AR Mixer has a driven mixer that can be removed quickly and easily. A major advantage when sterilizing the vessel in an autoclave or when sharing a single drive unit mechanism between vessels. A machine pad is welded into the bottom of the vessel (usually offset from a central bottom outlet valve). The impeller head is running by a drive unit mounts externally. Inside the impeller head there is a ring of special magnets encapsulated in a 316L stainless steel housing. Welded to this housing are the mixing blades which form an integral impeller head.

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Benefits

AR Inox has long recognised the benefits of Magnetic Mixer for applications in the pharmaceutical, biotechnological and cosmetic industries.

- Sanitary design to meet current pharmaceutical requirements.
- Mechanical seals are eliminated to give an aseptic design with no risk of contamination or leakage.
- Minimum working volume, which is not feasible with traditional mixers.
- Suitable for C.I.P. and S.I.P. in situ or autoclave.

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Performance

The geometry of the impeller allows an excellent axial pumping flow and very low shear rate for different applications:

- *mixing, diluting, maintaining in suspension, thermal exchange...*
- *Speed range 50 to 450 using an inverter. Depending on the size of vessel and impeller-head, levels of agitation can range from vigorous to gentle. Very small mixing volumes can be achieved.*
- *All product contact part are made in stainless steel AISI 316 L. The impeller head runs on a special grade tungsten carbide that has been manufactured under controlled conditions to achieve the correct size, grain structure and hardness. The running surface have a ground mirror finish.*
- *The blade shape and profile has been designed to optimise agitation levels and to direct a proportion of the cleaning fluid through the head.*

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Options

- A wide range of control can be supplied from a simple stand-alone panel to a fully integrated automation system.
- Magnetic proximity sensor can be fitted to measure the rotation of the impeller.
- Kit for vessel with jacket and insulation.
- Rotating blades directly welded to the magnetic head.
- Electro-polishing.

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Installation

The Magnetic Mixer is specially used in the pharmaceutical and biotechnological industries for ultra sterile application. The most obvious advantage of the magnetic mixer is that enables the total integrity of the tank.


The vessel shell is not penetrated and there are no mechanical seal. This also prevents toxic or high value media from leaking out of the vessel and thereby challenge the sterility of the vessel and its contents.

Suitable for viscosities up to 1000 cps the AR Magnetic Mixer offers great flexibility.

Bottom mounting also leaves the top of the vessel clear for easy vessel access.

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Mixing volume - Viscosity

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Key technical data

Model	Vigorous mixing	Gentle mixing
AR M4 DC 60W	10 LITRES	30 LITRES
AR M6 AC 0,18 KW	50 LITRES	150 LITRES
AR M8 AC 0,5 KW	150 LITRES	500 LITRES
AR M12 AC 0,75 KW	300 LITRES	2000 LITRES
AR M20 AC 1,5 KW	1000 LITRES	3000 LITRES
AR M32 AC 2,2 KW	2000 LITRES	6000 LITRES
AR M64D AC 5,5 KW	3000 LITRES	10000 LITRES